



MRI patient checklist

The system delivering HF10™ is approved for MRI scanning under certain conditions.

When a doctor requests an MRI scan:

1. *Tell the doctor who requested the MRI scan that you have an implanted spinal cord stimulation device.*
2. *Ask the doctor who implanted your system:*
 - *Can my system safely undergo the ordered MRI scan?*
3. *When scheduling your MRI appointment, provide the following information:*
 - *You have a Nevro spinal cord stimulator.*
 - *A detailed description of your devices, including the product number of the Implant, leads, and other implanted components.*
 - *The contact information of the doctor who implanted your device.*
 - *The MRI guidelines, which can be found at [HF10.com/resources/existing-patients](https://www.hf10.com/resources/existing-patients).*
4. *Take your patient ID card to the MRI appointment.*
5. *Prior to your MRI scan, make sure your stimulation is turned OFF. You can turn stimulation OFF with your Remote Control.*
6. *After the MRI scan, turn your system back ON with your Remote Control and contact your doctor to ensure that your HF10 is programmed to pre-MRI settings.*

Note: *Different models of implanted devices, including lead types, result in various approved scanning conditions. Please carefully review the MRI guidelines with your physician.*

