



# MRI patient checklist

*The system delivering HF10™ is approved for MRI scanning under certain conditions.*

## *When a doctor requests an MRI scan:*

- 1. Tell the doctor who requested the MRI scan that you have an implanted spinal cord stimulation device.*
- 2. Ask the doctor who implanted your system:*
  - Can my system safely undergo the ordered MRI scan?*
- 3. When scheduling your MRI appointment, provide the following information:*
  - You have a Nevro spinal cord stimulator.*
  - A detailed description of your devices, including the product number of the Implant, leads, and other implanted components.*
  - The contact information of the doctor who implanted your device.*
  - The MRI guidelines, which can be found at [www.Nevro.com/Manuals](http://www.Nevro.com/Manuals).*
- 4. Bring your patient ID card and Remote Control to the MRI appointment. If you don't have your patient ID card, please call your HF10 care team for assistance.*



*PTRC1000*



*PTRC2300*



*PTRC2500*

- 5. Prior to your MRI scan, make sure your stimulation is turned OFF. You can turn stimulation OFF with your Remote Control.*
- 6. After the MRI scan, turn your system back ON with your Remote Control and contact your doctor to ensure that your HF10 is programmed to pre-MRI settings.*

*Note: Different models of implanted devices, including lead types, result in various approved scanning conditions. Please carefully review the MRI guidelines with your physician.*

